


What to expect from your total knee replacement with patient-specific cutting guides



Patient guide



Just like your fingerprints, you are one of a kind.

Which is why every patient-specific cutting guide is designed to be just as unique as you are.

Total knee replacement

Years of research and development have gone into the creation of every element of the knee implant you are considering. Your implant's specialised shapes and the materials from which it is made have all been engineered, manufactured and tested to deliver the most accurate and reproducible outcomes possible.

That same level of care is also extended to the instruments that will be used during your surgery. After all, no matter how expertly designed an implant is, if it is not correctly placed in your body, it will never be able to perform at its peak. Consequently, neither will you.

Patient-specific cutting guides

It is important to note patient-specific cutting guides do not change the basic premise of a total knee replacement. To reconstruct the joint's functionality, your surgeon will remove the damaged bone and soft tissue in your knee and replace them with specially designed metal and plastic components. Patient-specific cutting guides allow your surgeon to select and position your new implant based on your anatomy before stepping into the operating room.

Your implant has been shaped to reflect the natural shape of a healthy knee. The back, or inner surfaces, of your implant consist of several flat surfaces which come together



Arthritic knee



Patient MRI



3D bone model
created from
patient's MRI

at precise angles. However, the process of shaping your bones to tightly fit these angles can be complex, invasive and time consuming because no two knee joints are exactly the same.

Patient-specific cutting guides provide a personalised solution to your surgeon to bridge the gap between the precisely engineered implant and your anatomy.

Advantages of patient-specific cutting guides

Shorter, less-invasive surgery


Patient-specific cutting guides reduce the number of steps required to perform total knee replacement surgery and eliminate the need to insert guide rods into the thigh bone. Both of which potentially shorten the amount of time a patient spends in the operating room and reduce some of the associated risks.¹⁻⁷

Faster, less-painful recovery

Along with a shorter, less-invasive procedure, other benefits may include the potential for a shorter hospital stay and less soft tissue disruption which can result in a faster road to recovery.

Implant longevity

Just like the tyres on your car, accurately aligning an implant can help improve its performance and extend its longevity.



Patient-specific cutting guide including surgeon-specific details



Surgeon requested bone cuts based on patient's anatomy



Properly placed implant

Knowing what to expect


Patient-specific cutting guides allow the surgeon to know what to expect with surgery. A preop plan is created for your case and identifies your implant size and placement that will be used. This plan is sent to the surgeon before the cutting guides are even made.

Why is the MRI so important?

- MRI is used to create a 3-D model of your anatomy
- The 3-D model allows the cutting guides to be designed to fit your specific bones
- The cutting guides based off the MRI identifies your implant size and positioning
- It is important no movement occurs during your MRI to ensure the highest quality and best fit of your cutting guides


The process

The inner shape of the cutting guide matches the outer shape of your bones. This hand-in-glove fit allows surgeons to make the precise bone cuts needed to position the knee implant in the optimal alignment.



Patient-specific cutting guide process

1. **Diagnostics:** Your surgeon will perform a diagnostic evaluation which will determine your need for a total knee replacement.
2. **MRI:** A short 8 – 10 minute MRI scan is performed several weeks prior to your surgery date.
3. **X-Ray:** A full leg standing X-Ray is also performed in conjunction with your MRI.
4. **Upload data:** Your MRI and X-Ray are sent to the manufacturer to begin the design process of your cutting guides.
5. **Pre-op planning:** A detailed plan of your surgery indicating your specific size is sent to your surgeon for review prior to your surgery date.
6. **Manufacturing:** Your cutting guides are manufactured and then shipped for your surgery.
7. **Surgery:** The cutting guides designed specifically for you are utilised for your surgery.
8. **Recovery:** After surgery, your surgeon will provide you with details about your recovery.



Disclaimer

Total knee replacement is not for everyone. Children, pregnant women, patients who have mental or neuromuscular disorders that do not allow control of the knee joint and morbidly obese patients should not undergo a knee replacement procedure. Total knee replacement is intended to relieve knee pain and improve knee function. However, implants may not produce the same feel or function as your original knee. There are potential risks with knee replacement surgery such as loosening, fracture, dislocation, wear and infection that may result in the need for additional surgery. Longevity of implants depends on many factors, such as types of activities and weight. This information, including postoperative care, is provided for educational purposes only. Smith & Nephew does not provide medical advice. In no event shall Smith & Nephew be liable for any damages whatsoever arising out of the use of or inability to use the expressed views.

1. Improved femoral component rotation in TKA using patient-specific instrumentation; Thomas J. Heyse, Carsten O. Tibesku; *The Knee*; November 2012.
2. Custom-fit total knee arthroplasty: our initial experience in 32 knees; Bali K, Walker P, Bruce W; *The Journal of Arthroplasty*, 2012; June 27 (6):1149-54.
3. The value of patient-matched instrumentation in total knee arthroplasty; Noble J, et al; *The Journal of Arthroplasty*; January 2012.
4. Frontal plane alignment after total knee arthroplasty using patient-specific instruments; Daniilidis K, Tibesku C.; *International Orthopaedics*; January 2013; 37 (1) 45-50.
5. 40420112 Smith & Nephew GENESIS® II DCF Surgical Technique.
6. Patient-Specific Versus Conventional Instrumentation for Total Knee Arthroplasty: Peri-Operative and Cost Differences; DeHaan AM, Huff, TW; *The Journal of Arthroplasty*; 2014.
7. Patient-specific instruments as a standard procedure in total knee arthroplasty: Logistics and postoperative radiological results in 70 patients, Drnek, D., Haffner, N., Sadjed A., Ritschl P.; *Case Reports in Clinical Medicine*; 2014.

Supporting healthcare professionals for over 150 years

Smith & Nephew Pty Ltd

85 Waterloo Road
North Ryde NSW 2113 Australia

Smith & Nephew Limited

Unit A, 36 Hillside Road
Wairau Valley Auckland 0627 New Zealand

°Trademark of Smith & Nephew

02986_AU V1 10/17 SN13612